



EYEWEAR OF YEAR 2016 AWARD WINNER (JAPAN)

OPTRA FAIR AWARD WINNER 2013

KIDS DESIGN AWARD 2014 (JAPAN)

GOOD DESIGN (JAPAN)

GOOD DESIGN (KOREA)

ODMA AWARDS 2017 (AUSTRALIA)

TOMATO GLASSES

BABY

KIDS

JUNIOR



NOV 2018

Our Story

Our journey started when I noticed that my 24 month old son's eyes were not aligning correctly. Our trip to the ophthalmologist confirmed Accommodative Esotropia.

His condition Hyperopia is a result of a small eye, which caused his eyes to strain so much that they crossed and eventually one eye became lazy.

The only correction was to wear thick convex lenses and to wear an eye patch to prevent a lazy eye known as Amblyopia.

The limited choice of frames on the market at the time with nose pads were only in metal. They were heavy and bulky and constantly kept falling off his face. Compliance was poor as my son refused to wear them due to the discomfort and weight. His condition worsened.

There was nothing available that could cater for a small child's head with a flat bridge. My disappointment was such that I was compelled to make his frames on my own.



This sparked the beginning of Tomato Glasses: a company and product born out of the affection and love I have for my child. His vision improved and with the help of surgery his vision is now normal.

I am glad I did not give up on him and ever so grateful that the Tomato Glasses helped him. I am also humbled to know that my passion and love for my son is now helping and supporting other children and parents across the globe that require custom fitted, lightweight and comfortable glasses. Children can now leave these spectacles on, with amazing success and great comfort. Tomato Glasses are frames that fit to a child's face rather than the child fitting the frame. We, Tomato Glasses, will continue to make children's frames so that we can make children and their parents happy.

THE
right
SIZE for
LITTLE
eyes




TOMATO GLASSES FEATURES



FIRM GLASSES

1. The frames are made using non-toxic materials that are approved by the FDA.
2. The frames are very flexible, display excellent elasticity and do not break easily.



COMFORTABLY FITTING GLASSES

1. The frames are incredibly light: only 6 to 9g in weight.
2. There is a cushion on the tip of the frames that matches the complex curves of the ears and will not create pressure points.

Weight by type

Baby 6.3g | Kids 7.0-9.5g | Junior 8.5-9.5g




NON-SLIP GLASSES

The bridge of the glasses matches the curve of the nose and therefore holds the glasses comfortably in place and prevents them from slipping down.




ADJUSTABLE TEMPLE LENGTHS

It is easy and simple to adjust the length of the temples. To shorten them, remove the ear tip and cut off the excess length. Replace ear tips and secure these in position using the screws. To lengthen the temples, pull the ear tips back and secure in place with the screws.



ADJUSTABLE NOSE PAD HEIGHTS

The nose pads can be attached to three points of the frame (five points on the baby range frames). As every child's face is unique, this feature allows the correct placement and height to be achieved for each child.



A STANDARD HEAD STRAP

All Tomato Glasses frames come with a head strap which attaches to the ear tips. They have been designed to prevent the glasses from slipping even with high powered lenses. The head strap provides security when running and jumping, making Tomato Glasses great for active children. They are also perfect for children with sensitive ears.





birth to
2yr s

Baby A

A very flexible frame for the utmost comfort which is highly adjustable to accommodate a baby's rapid growth.



The 35-13 frame is recommended for newborn babies up to approximately 6 months old while the 37-13 frame is for babies approximately 6 months to 1 year old and the 39-14 up to approximately 2 years old.

As a baby's face shape will grow during their early years, the TBAC range is highly adjustable and will accommodate these changes. The TBAC frame can be used with high powered lenses from -25 diopters through to +33 diopters. The baby frames can be adjusted along the horizontal and

vertical axis. With 5 positions available to attach the nose pads, the nose pad height can be adjusted to provide the correct position.

There are 3 different sized nose pads that can change the distance between the eyes and lenses. The flexible temples, comfortable ear tips and headband ensures the frame remains in place when babies are lying down or breast feeding.




Special features

POINT 1



There are 5 holes on the front frame to attach the nose pad.

POINT 2



There are 3 different sizes of nose pads that can decide the distance between the eyes and lenses.

POINT 3



Safe material, flexible and adjustable temples.

TBAC1 | Size 35□13 / 37□13 / 39□14




TBAC2 | Size 35□13 / 37□13 / 39□14



TBAC3 | Size 35□13 / 37□13 / 39□14



Kids Collections



These styles are designed to
colours with illustrations or s



Kids A

Oval frames

Kids B

Circular frames



d to appeal to young children using bright
r shapes on the sides of the temples.



2-8yrs



New

1-7yrs



Kids C

Rectangular frames

Kids D

Rectangular frames

KIDS A



TKAC1 | Size 37□14 / 40□14 / 43□15 / 45□16



TKAC25 | Size 37□14 / 40□14 / 43□15 / 45□16



TKAC3 | Size 37□14 / 40□14 / 43□15 / 45□16



TKAC26 | Size 37□14 / 40□14 / 43□15 / 45□16



TKAC4 | Size 37□14 / 40□14 / 43□15 / 45□16



TKAC29 | Size 37□14 / 40□14 / 43□15 / 45□16




TKAC7 | Size 37□14 / 40□14 / 43□15 / 45□16



TKAC31 | Size 37□14 / 40□14 / 43□15 / 45□16



TKAC8 | Size 37□14 / 40□14 / 43□15 / 45□16



TKAC11 | Size 37□14 / 40□14 / 43□15 / 45□16



TKAC19 | Size 37□14 / 40□14 / 43□15 / 45□16



TKAC12 | Size 37□14 / 40□14 / 43□15 / 45□16



TKAC27 | Size 37□14 / 40□14 / 43□15 / 45□16



TKAC13 | Size 37□14 / 40□14 / 43□15 / 45□16



TKAC28 | Size 37□14 / 40□14 / 43□15 / 45□16



TKAC14 | Size 37□14 / 40□14 / 43□15 / 45□16



TKAC32 | Size 37□14 / 40□14 / 43□15 / 45□16



TKAC15 | Size 37□14 / 40□14 / 43□15 / 45□16

KIDS B



TKBC1 | Size 40□14 / 42□15 / 44□16 / 46□17



TKBC7 | Size 40□14 / 42□15 / 44□16 / 46□17



TKBC2 | Size 40□14 / 42□15 / 44□16 / 46□17



TKBC8 | Size 40□14 / 42□15 / 44□16 / 46□17



TKBC3 | Size 40□14 / 42□15 / 44□16 / 46□17



TKBC9 | Size 40□14 / 42□15 / 44□16 / 46□17




TKBC4 | Size 40□14 / 42□15 / 44□16 / 46□17



TKBC21 | Size 40□14 / 42□15 / 44□16 / 46□17



TKBC5 | Size 40□14 / 42□15 / 44□16 / 46□17



TKBC11 | Size 40□14 / 42□15 / 44□16 / 46□17



TKBC17 | Size 40□14 / 42□15 / 44□16 / 46□17



TKBC12 | Size 40□14 / 42□15 / 44□16 / 46□17



TKBC18 | Size 40□14 / 42□15 / 44□16 / 46□17



TKBC13 | Size 40□14 / 42□15 / 44□16 / 46□17



TKBC19 | Size 40□14 / 42□15 / 44□16 / 46□17



TKBC14 | Size 40□14 / 42□15 / 44□16 / 46□17



TKBC22 | Size 40□14 / 42□15 / 44□16 / 46□17



TKBC15 | Size 40□14 / 42□15 / 44□16 / 46□17



TKCC1 | Size 40□14 / 42□15 / 44□16 / 46□17



TKCC6 | Size 40□14 / 42□15 / 44□16 / 46□17



TKCC2 | Size 40□14 / 42□15 / 44□16 / 46□17



TKCC7 | Size 40□14 / 42□15 / 44□16 / 46□17



TKCC3 | Size 40□14 / 42□15 / 44□16 / 46□17



TKCC8 | Size 40□14 / 42□15 / 44□16 / 46□17




TKCC4 | Size 40□14 / 42□15 / 44□16 / 46□17



TKCC9 | Size 40□14 / 42□15 / 44□16 / 46□17



TKCC5 | Size 40□14 / 42□15 / 44□16 / 46□17



TKCC11 | Size 40□14 / 42□15 / 44□16 / 46□17



TKCC16 | Size 40□14 / 42□15 / 44□16 / 46□17



TKCC12 | Size 40□14 / 42□15 / 44□16 / 46□17



TKCC17 | Size 40□14 / 42□15 / 44□16 / 46□17



TKCC13 | Size 40□14 / 42□15 / 44□16 / 46□17



TKCC18 | Size 40□14 / 42□15 / 44□16 / 46□17



TKCC14 | Size 40□14 / 42□15 / 44□16 / 46□17



TKCC19 | Size 40□14 / 42□15 / 44□16 / 46□17



TKCC15 | Size 40□14 / 42□15 / 44□16 / 46□17



TKDC1 | Size 38□17 / 40□17 / 42□17 / 44□17



TKDC5 | Size 38□17 / 40□17 / 42□17 / 44□17



TKDC2 | Size 38□17 / 40□17 / 42□17 / 44□17



TKDC6 | Size 38□17 / 40□17 / 42□17 / 44□17



TKDC3 | Size 38□17 / 40□17 / 42□17 / 44□17



TKDC7 | Size 38□17 / 40□17 / 42□17 / 44□17



TKDC4 | Size 38□17 / 40□17 / 42□17 / 44□17



TKDC8 | Size 38□17 / 40□17 / 42□17 / 44□17



TKDC11 | Size 38□17 / 40□17 / 42□17 / 44□17



TKDC15 | Size 38□17 / 40□17 / 42□17 / 44□17



TKDC12 | Size 38□17 / 40□17 / 42□17 / 44□17



TKDC16 | Size 38□17 / 40□17 / 42□17 / 44□17



TKDC13 | Size 38□17 / 40□17 / 42□17 / 44□17



TKDC17 | Size 38□17 / 40□17 / 42□17 / 44□17



TKDC14 | Size 38□17 / 40□17 / 42□17 / 44□17




TKDC18 | Size 38□17 / 40□17 / 42□17 / 44□17



Junior A

Sophisticated looking frame with a rectangular shape for children around 7 years old to adults looking for a more mature image.



TJAC1 | Size 47□17 / 49□17 / 51□17



TJAC6 | Size 47□17 / 49□17 / 51□17



TJAC2 | Size 47□17 / 49□17 / 51□17



TJAC11 | Size 47□17 / 49□17 / 51□17



TJAC3 | Size 47□17 / 49□17 / 51□17



TJAC12 | Size 47□17 / 49□17 / 51□17



TJAC4 | Size 47□17 / 49□17 / 51□17



TJAC13 | Size 47□17 / 49□17 / 51□17



TJAC5 | Size 47□17 / 49□17 / 51□17



7yrs to
adult

Junior B

Retro-style circular shaped frame projecting a gentle image for children aged 7 years old to adult.



TJBC4 | Size 46□17 / 48□17 / 50□18 / 52□18



TJBC12 | Size 46□17 / 48□17 / 50□18 / 52□18



TJBC5 | Size 46□17 / 48□17 / 50□18 / 52□18



TJBC13 | Size 46□17 / 48□17 / 50□18 / 52□18



TJBC6 | Size 46□17 / 48□17 / 50□18 / 52□18



TJBC14 | Size 46□17 / 48□17 / 50□18 / 52□18



TJBC7 | Size 46□17 / 48□17 / 50□18 / 52□18



TJBC15 | Size 46□17 / 48□17 / 50□18 / 52□18



TJBC9 | Size 46□17 / 48□17 / 50□18 / 52



Junior C

Rectangular shaped frame for a smart and intelligent image for children aged 7 years old to adult.





TJCC1 | Size 46□17 / 48□17 / 50□17 / 52□17 / 54□17.5



TJCC6 | Size 46□17 / 48□17 / 50□17 / 52□17 / 54□17.5



TJCC2 | Size 46□17 / 48□17 / 50□17 / 52□17 / 54□17.5



TJCC7 | Size 46□17 / 48□17 / 50□17 / 52□17 / 54□17.5



TJCC3 | Size 46□17 / 48□17 / 50□17 / 52□17 / 54□17.5



TJCC8 | Size 46□17 / 48□17 / 50□17 / 52□17 / 54□17.5



TJCC4 | Size 46□17 / 48□17 / 50□17 / 52□17 / 54□17.5



TJCC9 | Size 46□17 / 48□17 / 50□17 / 52□17 / 54□17.5



TJCC5 | Size 46□17 / 48□17 / 50□17 / 52□17 / 54□17.5



TJCC11 | Size 46□17 / 48□17 / 50□17 / 52□17 / 54□17.5



NOSE PAD FITTING KIT



Further research has indicated the need to cater for a variety of nose shapes. We have therefore developed a range of different sized nose pads, Baby fit (B1, B2, B3), Asian fit (A1, A2, A3) and Wide fit (W1, W2, W3). All nose pads are interchangeable on frames and the fitting kit stores all nine sizes of nose pad. This kit will allow optometrists and opticians to custom fit the frames to each child therefore offering the best customer service.

Our Reserch and Development team have created a Spacer Kit. It is designed to create space between the lenses and a child's eyelashes or face when using high powered lenses.




NOSE PAD FITTING KIT

SPACER KIT

DEFAULT NOSE PADS INSTALLED WITH TOMATO GLASSES FRAMES

size	35	37	38	39	40	42	43	44	45	46	47	48	49	50	51	52	54
TBAC	B2	B2		B2													
TKAC		A1			A1		W1		W2								
TKBC					A1	W1		W1		W2							
TKCC					A1	W1		W1		W2							
TKDC			A1		A1	W1		W1									
TJAC											W2		W3		W3		
TJBC										W2		W2		W3		W3	
TJCC										W2		W2		W3		W3	W3

THE RELATIONSHIP BETWEEN THE TYPES OF NOSE PADS AND FITTING RESULTS



If the frame sits too far away from a child's eyes, a wider nose pad should be used (e.g. replace a W1 for a W3 nose pad).
If the frame sits too close to eyelashes then the above spacer kit can be used to create a greater space.

FOR INFANTS (B1, B2, B3), PLEASE REPLACE IT WITH THE LOWER NUMBER.

TOMATO GLASSES CASSETTE AND ACCESSORIES

Accessories are included in the Tomato Glasses cassette. The Baby frames have more accessories in their pack.
Very convenient for storing parts.



FRAME PACK



BABY FRAME PACK

PRESIDENTS' MESSAGE



It has been many years since I made Tomato Glasses. Sometimes when I meet a child who wears Tomato Glasses on the street, I ask if the child feels any discomfort with the Tomato Glasses. We will always listen to the needs of children and parents and strive to make better glasses.


President Sungjoon Kim




CERTIFICATION



Durability test



Bridge test



Temperature stability test



Certificate of Europe



Registration



Korea Certificate

Our Awards


Tomato Glasses are an award winning specialized brand of children's glasses. We supply frames worldwide and we have a good reputation for design and quality from ophthalmologists and opticians around the world.



① EYEWEAR OF THE YEAR 2016 AWARD WINNER (JAPAN) / ② KIDS DESIGN AWARD 2014 / ③ Optra Awards 2013 1st / ④ Certificate of Good Design Selection / ⑤ ODMA AWARDS 2017 (AUSTRALIA) / ⑥ Toxin-Not Detected TEST REPORT(ear tips) / ⑦ Toxin-Not Detected TEST REPORT(nose pads) / ⑧ TEST REPORT(frame) Toxin-Not Detected / ⑨ Certificate Europe



⑩ No.10-1179157 Certificate of patent / ⑪ No.10-1186259 Certificate of patent / ⑫ Certificate of Trade Service Mark Registration / ⑬ Certificate of Trademark Registration(JAPAN) / ⑭ Certificate of Trademark Registration(USA) / ⑮ No.30-0652187 Certificate of patent / ⑯ No.30-0652189 Certificate of patent / ⑰ No.30-0652192 Certificate of Design / ⑱ Certificate Of Origin



TOMATO GLASSES

HEAD OFFICE

14, Gudeok-ro 232beon-gil, Seo-gu, Busan 49223, Korea
TEL: +82-51-244-8760 FAX: +82-51-244-8761 www.tomatoglasses.com

FACTORY

21, Haeyang-ro 57beon-gil, Youngdo-gu, Busan 49007, Republic of Korea
TEL: +82-70-7776-8760

JAPAN OFFICE

2-20 Hakatachuo Bld.2F, Nakagofukumachi, Hakata-ku, Fukuoka-shi, Fukuoka, 812-0035, Japan
TEL: +81-92-263-8166 FAX: +81-92-263-8167 www.tomatoglasses.co.jp

DISTRIBUTORS

Korea:	TOMATO GLASSES (+82) 51-244-8760
Japan:	TOMATO GLASSES (+81) 92-263-8166
United Kingdom:	TOMATO GLASSES UK (+44) 1395-271417
Europe:	TOMATO GLASSES EU (+44) 1395-271417
Australia:	TOMATO GLASSES AUSTRALIA (+61) 3-9448 8932
United States:	TOMATO GLASSES USA (+1) 718-309-0336
Canada:	TOMATO GLASSES CANADA (+1) 778-954-7880
Czech Republic:	TOMATO GLASSES CZECH REPUBLIC(+420) 605-747-864
Singapore:	TOMATO GLASSES (SINGAPORE) PTE LTD (+65) 6514 0730
Malaysia:	TOMATO GLASSES (SINGAPORE) PTE LTD (+65) 6514 0730
Philippines:	TOMATO GLASSES PHILIPPINES (+63) 920 932 6578
South Africa:	Danker Lenses (+27) 11-262-3980
Brazil:	Vis-à-vis Importação (+55) 11-26797902
New Zealand:	Obliquity limited trading (+64) 21-869-300
India:	Eleven2six Store (+91) 991-197-5008
Colombia:	Tienda Optica (+57) 4-313-8418
Taiwan:	Abby Optical Ltd (+886) 4-2471-1616 / GinDin International Co (+886) 2-2273-7785
Hong Kong:	Dreams Enterprise Ltd (+852) 3486-6882
Turkey:	MK OPTİK İÇ VE DIŞ TİCARET LİMİTED ŞİRKETİ (+90) 312 467 27 12
Indonesia:	PT.Optik Klinik Mata Nusantara (+62) 21-7581-6601